$$A = \pi (r_1 \times S_1 - r_2 \times S_2)$$

$$M = F \times r$$

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
Norton, M.  
An in vitro Evaluation of the Strength of an Internal Conical Interface Compared to a Butt Joint Interface in Implant Design  
Clin Oral Impl Res 1997; 8: 290-298

van Steenberghe, D., De Mars G., Quirynen M., Jacobs R., Naert I.  
A Prospective Split-mouth Comparative Study of Two Screw-shaped Self-tapping Pure Titanium Implant Systems  
Clin Oral Impl Res 2000; 11: 202-209

Hansson, S.  
The Implant Neck: Smooth or Provided with Retention Elements  
Clin Oral Impl Res 1999; 10: 394-405

Hansson, S., Werke, M.  
The Implant Thread as a Retention Element in Cortical Bone: The Effect of Thread Size and Thread Profile: A Finite Element Study  
J Biomechanics 2003; 36: 1247-1258

Abrahamsson, I., Berglundh, T.  
Tissue Characteristics at Microthreaded Implants: An Experimental Study in Dogs  
Clin Implant Dent Rel Res 2006; 8: 107-113

Wennström, J. et al.  
Implant-supported Single-tooth Restorations: A 5-Year Prospective Study  
J Clin Periodontol 2005; 32:567-574

Hansson S.  
Surface Roughness Parameters as Predictors of Anchorage Strength in Bone: A Critical Analysis  
J Biomechanics 2000; 33: 1297-1303

Hansson, S., Norton, M.  
The Relation Between Surface Roughness and Interfacial Shear Strength for Bone-anchored Implants. A Mathematical Model  
J Biomechanics 1999; 32: 829-836

Goffredsen, K., Nimb, L., Hjørting-Hansen E., Jensen, J.S., Holmén, A.  
Histomorphometric and Removal Torque Analysis for TiO<sub>2</sub>-Blasted Titanium Implants. An Experimental Study on Dogs  
Clin Oral Impl Res 1992; 3: 77-84

Ivanoff, C.J., Hallgren, C., Widmark, G., Sennerby, L., Wennerberg, A.  
Histologic Evaluation of the Bone Integration of TiO<sub>2</sub> Blasted and Turned Titanium Microimplants in Humans  
Clin Oral Impl Res 2001; 12: 128-134

Ellingsen, J. E., Johansson, C., Wennerberg, A., Holmén, A.  
Improved Retention and Bone-to-Implant Contact with Fluoride-Modified Titanium Implants  
Int J Oral Maxillofac Implants 2004; 19: 659-666

Cooper, L., Zhou, Y., Takebe, J., Guo, J., Abron, A., Holmén, A., Ellingsen, J.E.  
Fluoride Modification Effects on Osteoblast Behavior and Bone Formation at TiO<sub>2</sub> Grit-blasted c.p Titanium Endosseous Implants  
Biomaterials 2006; 27: 926 – 936

Isa Z.M., Schneider, G., Zaharias, R., Seabold, D., Stanford, C.M.  
Effects of Fluoride-Modified Titanium Surfaces on Osteoblast Proliferation and Gene Expression  
Int J Oral Maxillofac Implants 2006; 21: 203-211

Stanford, C., Johnson, G., Fakhry, A., Gartton, D., Mellonig, J., Wagner, W.  
Outcomes of a Fluoride Modified Dental Implant One Year after Loading in the Posterior-Maxilla when Placed with the Osteotome Surgical Technique  
Applied Osseointegr. Res. 2006; 5: 50 – 55

Schliephake, H., Hüls A., Müller, M.  
Early Loading of Surface Modified Titanium Implants in the Posterior Mandible – Preliminary Results  
Applied Osseointegr. Res. 2006; 5: 56 – 58

Palmer, R.M., Palmer, P. J., Smith, B.J.  
A 5-Year Prospective Study of Astra Tech Single Tooth Implant  
Clin Oral Impl Res 2000; 11: 179-182

Puchades-Roman, L., Palmer, R.M., Palmer, P.J., Howe, L.C., Ide M., Wilson, R.F.  
A Clinical, Radiographic, and Microbiologic Comparison of Astra Tech and Brånemark Single Tooth Implants  
Clin Impl Dent Rel Res 2000; 2: 78-84

Abrahamsson, I., Berglundh, T., Wennström, J., Lindhe, J.  
The Peri-implant Hard and Soft Tissues at Different Implant Systems  
Clin Oral Impl Res 1996; 7: 212-219

Wennström, J., Ekestubbe, A., Gröndahl, K., Karlsson, S., Lindhe, J.  
Implant-supported Single-tooth Restorations: A 5-Year Prospective Study  
J Clin Periodontol 2005; 32:567-574

De Kok, I.J., Chang, S.S, Moriarty, J.D., Cooper, L.F.  
A Retrospective Analysis of Peri-Implant Tissue Responses at Immediate Load/ Provisionalized Microthreaded Implants  
JOMI 2006, 21:405-412

Norton, M.  
Multiple Single-Tooth Implant Restorations in the Posterior Jaws: Maintenance of Marginal Bone Levels with Reference to the Implant-Abutment Microgap  
Int J Oral Maxillofac Implants; in press.

Abrahamson, I., Berglundh, T., Lindhe, J.  
Soft Tissue Response to Plaque Formation at Different Implant Systems. A comparative study in the Dog.  
Clin Oral Impl Res 1998;9: 73-79